

Grid photopatterning: customized protocol

Responsible person: Raaja

Autors/Illustrations: Carlotta, Šárka

A. Reagents and Equipment

Equipment

- PRIMO® setup
- TEM grid
- Parafilm
- Glass slides
- Wet chamber
- Plasma cleaner or glow discharge system (optional depending on the passivation)
- Optional: hot plate, Al-foil




Summary

- A. Reagents and Equipment
- B. Patterning using PEG-SVA passivation and PLPP photoinitiator
- D. Coating with protein of interest
- B., D. Illustrated protocol
- C. Primo2
- D. Customized shapes/templates for Leonardo

Reagents

- Protein solution of interest
- 1x PBS, mQ H₂O, MeOH/EtOH
- HEPES pH 8.3 – 8.5 at 0.1 M ... **! critical:** needs to be check at each exp. **!**
- mPEG-SVA (mPEG-Succinimidyl Valerate, MW 5 000, Laysan bio)
- Poly-L-Lysine (PLL, Merck/ Sigma-Aldrich, ref: P8920)
- oriPLPP ... PLPP Gel (Alvéole)
- oursPLPP ... 14.5 mg/mL in MeOH
- mPEG-SVA is a PEG (5kDa) coupled with a Succinimidyl Valerate (SVA) ester. This ester is used for PEGylation of primary amine of the polylysine

HEPES 0.1M; pH 8.3 – 8.5 :
preferably freshly prepared
0.59 g ... 25 mL mQ H₂O
1.19 g ... 50 mL mQ H₂O
2.38 g ... 100 mL mQ H₂O
Adjust pH with NaOH




B. Patterning using PLPP Gel photoinitiator and PEG-SVA passivation

before you start – check Leonardo software, if calibration of Primo is alright (viz. Note)

B1. Passivation protocol with PEG-SVA


- Check pH of HEPES (freshly prepared or only a few days in advance)
- Glow discharge TEM grid
- Place grid onto parafilm in wet chamber (carbon side top)
- Add PLL: 100 µg/1 mL of HEPES (10 µL of stock + 90 µL HEPES) , 20 - 30 µl/grid, 30 min
- **!** Rinse 3 x 5 min with large volume of HEPES aprox. 200 - 500 µL (pH 8.3 – 8.5)
– pipet up and down several times during each wash
- Prepare a fresh solution of PEG-SVA 100 µg/mL (70 µg/mL - 100 mg/mL) in HEPES; during the last rinsing
- **Note:** As the half-life of the SVA ester is 10 minutes at pH 8.5, the solution must be prepared just before use.
. If the passivation is not strong enough, you can increase the concentration up to 100 mg/mL.
- **Recommendation:** Change parafilm (support, I do not do that and works well)
- Add PEG-SVA: 200 - 500 µl/grid, 60 min
- Rinse with 3 x 5 min mQ H₂O; pipet up and down several times during each wash



Stopping point: The passivated grid can be kept dried for few weeks at 4°C. Clean and sealed box on parafilm.

B2. Photopatterning with PRIMO using PLPP photoinitiator

- (**OriPLPP**: Prepare a solution of PLPP gel in mQ H₂O:
1 µl of PLPP gel / 2 µl of mQ H₂O;)
*Wait for the complete drying of the solution, a transparent gel will form
- **OursPLPP**: aliquot from Rajaa; add 66 µL of MeOH or EtOH; 3 µl
- Place the grid upside down on the glass slide (gold side top)
- Add a solution of PLPP : 3 µl/grid before patterning



Tip: Drying steps
Glass slide with grids can be placed on heating plate on Al-foil: 70°C; approx. 30 - 45 min

Note: At this point – in the meantime, you can start Primo2 and Leonardo software and proceed calibration on another glass slide (**C1 – C3; detailed protocol page 4**). Briefly:

- Focus on the surface of the glass slide
- Load and lock your pattern on Leonardo software
- Detection functionality or the TEM grid functionality (Check tutorial video at <https://www.alveolelab.com/tutorials/micropatterning-tem-grid/>)
- The typical UV doses : Primo2 (LED, 365nm): 20mJ/mm⁻². It might depend on the density of PLPP Gel, substrate, etc ... dose adjustment might be needed.
- Launch patterning sequence – then continue with your grids



Tip: Usefull numbers in Leonardo 1
Pixel size of pattern: 200 x 200 px
Laser: 25%
Grid diameter: 3000 µm
View: 36 square (200 mesh)


D. After Patterning **Note:** do not stop, finish until cell seeding

- Rehydrate the slide with PBS for another 5 – 10 minutes (30 min).
- **!** Rinse profusely with large volume PBS: 3 x 5min

D1. Coating with protein of interest

- Coat with protein of interest according to its specific procedure/your experience:
e.g. Incubation time 5 – 60 (120 min), 20 – 100 mg/mL
- **!** Rinse profusely with PBS without drying:
 - 3 x 5 min, large volume
 - keep PBS on grids for easier handling, ideally for less than 10 min

D1. Seed cells



- Glow discharge TEM grid
- Place grid on to parafilm in wet chamber

HEPES pH 8.3 – 8.5

- Prepare a solution of PLL in HEPES: 100 µg/mL (10 µL of PLL + 90 µL of HEPES)
- Add PLL: 20 - 30 µl/grid; 30 min
- Rinse 3 x with HEPES (200 - 500 µL)
- Prepare a solution of PEG-SVA in HEPES: 100 µg/mL (options: 70 µg/mL – 100 mg/mL)
- Add of PEG-SVA: 200 (500) µl/grid, 60 min
- Rinse 3 x with x mQ H₂O;
- The last drying before PLPP: quickly and very gently dip grid on Kimtech

Stopping point: The passivated grid can be kept dried for few weeks at 4°C. Before PLPP: Rehydrate the slide with PBS

- Place grid upside down on a glass slide oriPLPP
- Prepare solution oriPLPP 1 µl of PLPP gel / 2 µl of mQ H₂O;
- Wait for the complete drying of the solution (70°C, 1h, in dark)
- oursPLPP:
 - aliquot from Rajaa in MeOH/EtOH, directly use 3 µl/grid before pattern, w/o drying step

- C. Primo2:** pages 4 – 5
- Pattern your grids

- Rehydrate the slide with PBS (mQ H₂O) for 5 – 10 min directly on the glass slide (I usually fill glass slide full of PBS and leave it for 30 min)
- After rehydration: you can gently transfer grids somewhere ...
- Rinse profusely with PBS: 2 x 5 min

(I do that on Parafilm – on bench; can be done in the cell culture hood)

- Coat with protein of interest according to its specific procedure/your experience
- Rinse profusely with PBS without drying
Note - washing: 10-cm dish with parafilm: place grids onto parafilm, pour approx.. 10 mL PBS and pipet with pipetboy up and down several times, leave it under PBS for 5 min
do it twice

Timing for 12 grids

Leonardo cal. : 30 min
Glow discharge: 30 min
PPL: 30 min
Wash: 20 min
PEG: 60 min
Wash: 20 – 30 min
Primo + PLPP: 10 min/grid
Wash: 45 min
Protein of interest: 30 – 120 min
Wash: 30 min
Seeding: 60 min

- (sterilization of grids: UV)
- Seed cells

C. Primo2

C1.

- Open Micromanager
- Select: PRIMO
- Upper bar → Leonardo → Leonardo

C2. Controlle microscope:

- Magnification for grids: 20x
- EP1: mirror 1
- EPI2: empty/none
- DIA : LED intensity 50
- DIA : shutter OFF
- ZDC: in







C3. Leonardo software


I. Calibration

- Preparation
 - Magnification: 19.793
 - Laser power: 3
 - Lens 450
- Alignment
 - Done by Rajaa
- Focus
 - Start operation
 - Exposure: 25 ms
 - Laser: 25%
 - 'PRIMO' logo: focus aprox. 5000
focus on the smallest feature (ctrl + mouse wheel)
 - ZDC: in
 - DIA : shutter OFF







Tip: If the calibration was not done during PLPP, it is possible to draw highlighter line next to grids; but **not recommended**; the highlighter line can be smudged by PBS





II. Patterning

- I. Control panel: DIA: shutter ON
- II. Leonardo: upper bar: Laser 100%
- a. Find centre of grid – place it in the middle of view (ruler tool)
- b. Focus on squares/holes
- c. Set bright light at correct intensity (the min/max of the histogram) - to have a “black & white” image: meshes are completely white and their contours black.
- d. Wizard → TEM grids → NEXT
- e. Select a shape
- f. * Diameter of ROI, (dose), shape size
→ add grid → ROTATE → ☒ keep ratio → LOCK
- g. FINISH
- h. Wait for ►
- i. Press ►
- j. Wait until finish
- k. Change grid → repeat: a. → j.




- Orange dots: found after (g) FINISH
- ►
- Blue dots: in progress
- Yellow dots: done

- *diameter of ROI : variable,
grid ... 3000 μm
the smallest area to be pattern = field of view ... square; 750 μm
→ for smaller ROI – it is possible to change
the ratio (to have rectangle ...)

→ LOCK

dose: during calibration



D. Customized shapes/templates for Leonard

I. Draw shape(s) in suitable software

- e.g. Illustrator, Inkspace, Affinity
Export in as:
 - Black shape
 - White background
 - TIFF
- Fiji ImageJ
 - Import to Fiji – drag TIFF
 - Edit (upper bar) → Invert
 - Image (upper bar) → Type → 8-bit
 - Image (upper bar) → Scale → Width 200px; Height 200px

II. Transfer to microscope computer

- via your Citrix there: copy/paste
- Note: it doesn't work from Drop-in and there is no internet connection

III. Transfer to Leonardo

- Put it into:
Micromanager 1.4 'folder' → Patterns

SaS – calculations

12 grids → 200 µL/each → 2400 µL

PPL:

12 grids: 40 µL + 360 µL

33 µL/grid

PEG:

12 grids: 2.4 mg + 2.40 mL

200 µL/grid

PLPP:

Ours from Rajaa: add 3 µL immediately before each grid

can be add also before start and just continue

Both ways works – there is only one problem – tilting of grids – saturation can help